



## Glasswall TQ

Durable, Highly  
Abrasion and Chemical  
Resistant Epoxy  
Coating

## SeamTek™ Glasswall TQ Epoxy Wall Coating System



## SeamTek™ Glasswall TQ Features

- Extremely Durable
- No amine blush–no frosting
- Color Stable
- Decorative Wall Surface
- Superior Impact Resistance
- Excellent adhesion to concrete
- Good Workability – easy to spread
- 100% solids – Solvent Free
- Low VOC -Low Odor
- High Taber Resistance
- Low Flammability

## Life Science Products

124 Speer Road, Chestertown, MD 21620  
[www.lspinc.com](http://www.lspinc.com) | 800-638-9874 | [info@lspinc.com](mailto:info@lspinc.com)

© 2022 Life Science Products, Inc.

# SeamTek™ Glasswall TQ

## General Description

SeamTek™ GlassWall TQ is a three-component 100% solids, low-odor, low-VOC epoxy resin wall covering system. It chemically cures to form a highly abrasion-resistant wall system when high-performance interior surface finishes are required. SeamTek™ Glasswall TQ is designed to exhibit decorative characteristics, impact resistance, chemical resistance, and durability.

It has been formulated as a mortar to be a nominal 1/8" thick. This is the same resin and aggregate formulation as an epoxy quartz floor. The basic resin component (SR101) is formulated to work with several selected hardeners to serve as a highly durable matrix. The seal coat hardener (CRH405) assures outstanding chemical resistance when compared to other CR epoxy products.

SeamTek™ GlassWall TQ can be installed over masonry, poured concrete and hardboard surfaces to provide a seamless aesthetic surface. Unless otherwise specified, SeamTek™ GlassWall TQ will follow the contour(s) of the existing substrate and cannot be used as a stand-alone system to correct such problems. SeamTek™ GlassWall TQ is designed to be 1/8" thick, which is referred to as a mortar mix system.

## SeamTek™ Glasswall TQ Details and Properties

**Specific Gravity** - ASTM D 70, Fisher #3-247

Pycnometer – 1.25

**Weight:** ASTM E 201 – 10.7 lbs./gal

**Non-Volatile Content:** ASTM D 1353, 18 hrs. at 200°F (93°C) – 100%

**Viscosity, cps** - ASTM D 1475 77°F (25°C) 8,850 cps

**Flash Point TCC Minimum** - Seta Flash – Greater Than 2000 F (930 C)

**Solvent Odor** – ASTM D1296 – Extremely Low

**Pot Life** – 20 Minutes at 720 F (220 C) & 50% R.H.

**Flammability** - ASTM D 635 - Self-Extinguishing

## Mixing:

1. Apply 101 epoxy primer to the substrate and seed with fine aggregate. Allow to cure.
2. Basic mortar mix: SR101 and SH101 in 2:1 mix ratio as above. Mix 1 ½ quarts (three pints) of mixed liquid with 11 quarts of color quartz or 50 mesh sand. Trowel on wall @ 12 – 13 square feet per gallon of mortar mix.
3. Grout Coat: Mix neat resin SR101 and hardener CRH405 in a 2:1 ratio and apply at a spread rate sufficient to fill the surface porosity. Repeat in successive applications until "holidays" disappear. This may require successive days.
4. Seal Coat: Mix neat resin SR101 and hardener CRH405 in a 2:1 ratio and apply at a spread rate sufficient to leave an even sealed surface.

## Storage and Handling:

Because SeamTek® GlassWall TQ has a flash point above 200°F (93°C), transportation, storage and handling are less restricted.

## Limitations:

SeamTek® GlassWall TQ must not be used to bridge moving cracks or joints. Non-moving cracks or joints that must be over coated require rigid repairs. Surface or air temperature must be between 65°F minimum and 75°F maximum and relative humidity below 80%. Lower temperatures will extend cure time and higher temperatures will reduce pot and work life.

## Product Health and Safety Information:

Refer to container labels and Material Safety Data Sheets available from LSP for health, safety and environmental information.

*If necessary, call LSP at (800) 638-9874.*

**Life Science Products have been in demand by these and other highly respected institutions:**

Bristol Meyer Squib | Children's Mercy | Cleveland Clinic | CalTech Univ. | Dana Farber | Duke University  
Emory University | F.D.A. | Harvard University | M.D. Anderson | NIH | Novartis | Northwestern University  
Ohio State U. | Pfizer | Princeton University | Regeneron | University of North Carolina | Yale University